

GAMUT

DIN 02241575

Product	Multi-purpose disinfectant, decontaminant and cleaner.
Description	A ready-to-use, multi-purpose, broad-spectrum safe and compatible solution for use on environmental surfaces and inanimate objects.
Function	GAMUT will effectively clean and decontaminate critical and semi-critical instrumentation prior to appropriate sterilization or high-level disinfection. The solution can also be used to disinfect non-critical devices such as stainless steel or plastic trays.
Packaging	4 L economy size with complimentary empty spray bottle.

Features

Benefits

Ready-to-use	GAMUT is supplied ready-to-use in an economical 4 L bottle. No dilution is required. Stable 3-year shelf life affirms efficacy.
Multi-Purpose	Non-toxic GAMUT can be used as a hard surface disinfectant, instrument decontaminant, ultrasonic or presoak solution.
Rapid & Effective	The synergistic combination of alcohol and quaternary ammonium means the solution is bactericidal, virucidal and fungicidal in just 1 minute. Complete Tuberculocidal disinfection is achieved in 10 minutes.
Broad Spectrum Anti-Microbial Agent	Maximum disinfection against most disease and odour causing organisms: Staphylococcus aureus Pseudomonas aeruginosa Herpes simplex virus type blend 2 Salmonella choleraesuis Trichophyton mentagrophytes Polio virus type 1 and 2 Coxsackie virus Mycobacterium Tuberculosis HIV-1 (AIDS Virus) MRSA VRE
Safe to Use	Biodegradable GAMUT is safe for use on tables, lights, plastics, stainless steel, Plexiglas, glass, rubber, and other non-porous surfaces. GAMUT is especially useful in eliminating environmental cross-contamination from mobile medical equipment.
Controls Odours	Destroys mildew, mold and fungi.
Fine Aerosol Mist	Effective penetration of fibrous materials, cracks and crevices.
DIN Registered	Assurance that all claims stated on label have been verified.